



### POSITIONING PRODUCTS RANGE



# **STRAIGHT UP** THE BEST QUALITY

KINETIC BALANCE'S MISSION IS TO DEVELOP AND PRODUCE HIGH QUALITY REHABILITATION AIDS THAT ARE DISTINCTIVE IN APPEARANCE AND FUNCTIONALITY.

CHEST POSITIONING	6/7
CHEST POSITIONING FLEX CONTROL	8/9
PRODUCT CHARACTERISTICS	10/11
2-POINT HIP POSITIONING	12
4-POINT HIP POSITIONING	13
PELVIC POSITIONING	- 14
ANKLE POSITIONING	15
ABDUCTION POSITIONING	16
NEOPRENE CHEST BELT X	18
NEOPRENE MULTI-PURPOSE BAND	19
SHOWER BELT	19
PRODUCT OVERVIEW & ITEM NUMBERS	20/21

KINE COMME

# **PRODUCT CHARACTERISTICS**

### UNIQUE FEATURES BY KINETIC BALANCE

'Kinetic Balance' has been providing its fresh take on positioning since its inception in 2004. Kinetic Balance products are detail-oriented and set themselves apart from the competition thanks to their unique features. Kinetic Balance sets new standards time and again. First with the introduction of waterresistant materials and covered fasteners, with our flex-control vests and continuous improvements to both the construction, look and feel.



WATER AND DIRT REPELLENT



TRIMSTRIP



KB'S ORIGINAL, COVERED FASTENERS



QUALITY CONTROL



GRIP PRINT



ROUND WEBBING ENDS





SUPER-DURABLE STITCHING



PREMIUM QUALITY MATERIALS



CLEAR SIZE INDICATION



ALL FASTENERS INCLUDED

### AND MORE ...

- SOFT NEOPRENE EDGE
- HIGH QUALITY LIMESTONE NEOPRENE
- WRINKLE FREE WEBBING
- ISO 9001:2015 PRODUCTION FACILITY



ACCORDING TO DIRECTIVE 93/42/EEC



MATERIAL REINFORCEMENTS

# ORIGINAL KINETIC BALANCE IDEAS FOR POSITIVE POSITIONING.


# CHEST POSITIONING

ALWAYS USE IN COMBINATION WITH A HIP BELT

### WOMEN CHEST POSITIONING

ARTICLE NUMBER: 5.71014.X SIZE: S & M





### MEN CHEST POSITIONING

ARTICLE NUMBER: 5.71012.X SIZE: XS TO L



### **CHEST POSITIONING XS**

ARTICLE NUMBER: 5.71012.0 SIZE: XS



### MEN CHEST POSITIONING ZIPPER

ARTICLE NUMBER: Z5.71022.X SIZE: S & M





	A	В	С	D
XS	190	210	140	600
S	220	300	210	600
JUNIOR	220	300	220	600
М	250	350	230	600
L	280	400	250	600
		All 10		nto in mm

# CHEST POSITIONING FLEX-CONTROL

ALWAYS USE IN COMBINATION WITH A HIP BELT





THE FLEX-CONTROL LETS YOU DETERMINE THE LEVEL OF FLEXIBILITY TO SUIT YOUR NEEDS THE TRIM STRIP ON THE ADJUSTABLE STRAP MAKES CORRECT ATTACHMENT EVEN SIMPLER


### WOMEN CHEST POSITIONING FLEX-CONTROL

ARTICLE NUMBER: 5.71034.X SIZE: S & M

#### MEN CHEST POSITIONING FLEX-CONTROL

ARTICLE NUMBER: 5.71032.X SIZE: S TO XL





М	400	50	100	100		
L	480	50	120	120		
XL	540	70	150	140		
		All measurements in mm				

#### **H-CHEST POSITIONING FLEX-CONTROL**


ARTICLE NUMBER: 5.81001.X SIZE: S TO XL



# MEASURING

### **CHEST POSITIONING**





TRUNK LENGTH MEASURED, WHEN SEATED, FROM INCISURA JUGULARIS TO SIAS

### BELTS, ABDUCTION AND PELVIC POSITIONING



HIP MEASURED FROM SIAS TO SIAS

### **ANKLE POSITIONING**



CIRCUMFERENCE MEASURED OVER MALLEOLI

# **RELEASE TYPES**

SIDE RELEASE

SIDE RELEASE

### **CENTER RELEASE**

**CENTER RELEASE** 

### **PIN RELEASE**

### **AIRPLANE BUCKLE**



# **2-POINT HIP BELT**

#### HIP BELT DOUBLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70001.X



#### HIP BELT CENTER RELEASE SINGLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70003.X



### HIP BELT PIN RELEASE DOUBLE ADJUSTABLE

MAAT: XS TO XL ARTICLE NUMBER: 6.70005.X



#### HIP BELT SINGLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70002.X



#### HIP BELT CENTER RELEASE DOUBLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70004.X



#### HIP BELT AVIATION RELEASE SINGLE ADJUSTABLE MAAT: S TO XL

ARTICLE NUMBER: 6.70006.X



# **4-POINT HIP BELT**

#### 4-POINT HIP BELT SINGLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70042.X



#### 4-POINT HIP BELT CENTER RELEASE DOUBLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70043.X

#### 4-POINT HIP BELT DOUBLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70044.X



# 4-POINT HIP BELT PIN RELEASE DOUBLE ADJUSTABLE

SIZE: XS TO XL ARTICLE NUMBER: 6.70045.X






### 4-POINT HIP BELT CENTER RELEASE ADJUSTABLE

MAAT: XS TO XL ARTICLE NUMBER: 6.70047.X





	A	В	С	D	E
XS	50	175	40	600	500
S	70	175	50	650	500
м	75	220	55	750	600
L	80	280	55	800	600
XL	95	330	60	800	700

# **PELVIC POSITIONING**

#### PELVIC POSITIONING NON-STRETCH

ARTICLE NUMBER: 7.70001.X SIZE: S TO XL

10



ARTICLE NUMBER: 7.70002.X SIZE: S TO XL





	A	В	С	D	E	F
S	180	180	45	45	550	400
М	220	220	50	50	600	450
L	300	300	60	60	650	500
XL	400	400	80	80	750	550

# **ANKLE POSITIONING**

#### ANKLE POSITIONING QUICK RELEASE CLOSURE

ARTICLE NUMBER: 8.70002.X SIZE: XS TO L




ARTICLE NUMBER: 8.70001.X SIZE: XS TO L

#### ALWAYS USE IN COMBINATION WITH A HIP BELT



Contour of the malleolus lateralis/medialis

O < 18cm = XS
< 22cm = S
< 26cm = M
> 26cm = L




	A	В	С	D	E	F	G
XS	80	65	55	45	80	210	210
S	110	75	60	55	110	210	210
м	140	95	65	60	140	260	260
L	170	110	65	60	170	260	260

# ABDUCTION POSITIONING

BALA



ARTICLE NUMBER: 9.70001.X SIZE: XS TO L



# ISO 9001:2015 MANUFACTURING STANDARD

ISO 9001 IS THE INTERNATIONAL STANDARD THAT SPECIFIES REQUIREMENTS FOR A QUALITY MANAGEMENT SYSTEM (QMS). ORGANIZATIONS USE THE STANDARD TO DEMONSTRATE THE ABILITY TO CONSISTENTLY PROVIDE PRODUCTS AND SERVICES THAT MEET CUSTOMER AND REGULATORY REQUIREMENTS.

OUR MANUFACTURING FACILITY IS COMPLIANT WITH THIS STANDARD. THIS ASSURES THE SUPERIOR QUALITY OF ALL OUR PRODUCTS, SO WE CAN AND WILL CONTINUE TO PROVIDE YOU WITH THE BEST POSITIONING PRODUCTS AVAILABLE.

NOTES:

# ASSORTED POSITIONING PRODUCTS

### **NEOPRENE CHEST BELT X**

THE NEOPRENE CHEST BELT X IS A BIG STEP FORWARD IN SOFT COMFORT CHEST POSITIONING

ARTICLE NUMBER: 402X

SIZE: XS TO L

	A	В	С	
XS	310	80	500	
S	400	100	500	
М	500	125	500	
L	570	130	500	
	All measurements in mm			



### NEOPRENE MULTI-PURPOSE BAND



### **SHOWER BELT**

THE SHOWER BELT IS MADE OF THE BEST MATERIALS AVAILABLE. IT ENSURES COMFORT AND FIT WHEN WET







# **PRODUCT OVERVIEW & ITEM NUMBERS**

### **CHEST**

Men's chest positioning				
5.71012.0	XS			
5.71012.1	S			
5.71011.0	Junior			
5.71012.2	М			
5.71012.3	L			

Woman's chest positioning				
5.71014.1 S				
5.71014.2	М			

Chest positioning with zipper				
Z5.71022.1	S			
Z5.71022.2	М			

Men's chest positioning, flex-control					
5.71032.1	S				
5.71032.2	М				
5.71032.3	L				
5.71032.4	XL				

Woman's chest positioning, flex-control		
5.71034.1	S	
5.71034.2	М	

H-positioning, flex-control, non stretch		
5.81001.1	S	
5.81001.2	М	
5.81001.3	L	
5.81001.4	XL	

Neoprene Chest Belt X		
4020	XS	
4021	S	
4022	М	
4023	L	

### **PELVIC**

Pelvic positioning non-stretch		
7.70001.1	S	
7.70001.2	М	
7.70001.3	L	
7.70001.4	XL	

Pelvic positioning stretch		
7.70002.1	S	
7.70002.2	М	
7.70002.3	L	
7.70002.4	XL	

## **ABDUCTION**

Abduction positioning		
9.70001.0	XS	
9.70001.1	S	
9.70001.2	М	
9.70001.3	L	

## 

Ankle positioning velcro closure		
8.70001.0	XS	
8.70001.1	S	
8.70001.2	М	
8.70001.3	L	

Ankle positioning quick release closure		
8.70002.0	XS	
8.70002.1	S	
8.70002.2	М	
8.70002.3	L	

### MULTI-PURPOSE

Neoprene multi-purpose band		
5002	M	

### HIP

2-point hip belt double adjustable		
6.70001.0	XS	
6.70001.1	S	
6.70001.2	М	
6.70001.3	L	
6.70001.4	XL	

2-point hip belt single adjustable		
6.70002.0	XS	
6.70002.1	S	
6.70002.2	М	
6.70002.3	L	
6.70002.4	XL	

2-point hip belt centre release single adjustable		
6.70003.0	XS	
6.70003.1	S	
6.70003.2	М	
6.70003.3	L	
6.70003.4	XL	

2-point hip belt centre release double adjustable			
6.70004.0	XS		
6.70004.1	S		
6.70004.2	М		
6.70004.3	L		
6.70004.4	XL		

2-point hip belt pin release double adjustable			
6.70005.0	XS		
6.70005.1	S		
6.70005.2	М		
6.70005.3	L		
6.70005.4	XL		

2-point hip belt aviation release single adjustable		
6.70006.1	S	
6.70006.2	М	
6.70006.3	L	
6.70006.4	XL	

4-point hip belt single adjustable		
6.70042.0	XS	
6.70042.1	S	
6.70042.2	М	
6.70042.3	L	
6.70042.4	XL	

4-point hip belt centre release		
6.70043.0	XS	
6.70043.1	S	
6.70043.2	М	
6.70043.3	L	
6.70043.4	XL	

4-point hip belt double adjustable		
6.70044.0	XS	
6.70044.1	S	
6.70044.2	М	
6.70044.3	L	
6.70044.4	XL	

4-point hip belt pin release		
6.70045.0	XS	
6.70045.1	S	
6.70045.2	М	
6.70045.3	L	
6.70045.4	XL	

4-point hip belt centre release adjustable		
6.70047.0	XS	
6.70047.1	S	
6.70047.2	Μ	
6.70047.3	L	
6.70047.4	XL	

## SHOWER

Shower belt		
3001	S	
3002	М	
3003	L	



# CLOTHING & BAGS FOR SEATED PEOPLE



WWW.KINETIC-BALANCE.COM





POSITIONING PRODUCTS BV HELLINGWEG 11-C 2583 DZ, THE HAGUE THE NETHERLANDS T: 0031 (0)70 217 09 08 E: INFO@KINETIC-BALANCE.COM E: INFO@KINETICBALANCE.NL



WWW.KINETIC-BALANCE.COM